



welcome to sheffmed

We are proud to be different. Our customer driven approach together with a team orientated work ethos, means that you receive a friendly & reliable service. We are an ISO registered company and all our products carry the CE mark. Our team are committed to helping you to provide the best quality care for your patients, particularly with the new budget restraints and of course the ever changing environment. Please do not hesitate to contact us if you require any additional information.



Canine Operatives:
Flyte, Dotty, Harry & Wispa



Office Based:
Sonia, Liam & Maxine



Sales Team:
Ian, Phill, Philip & Alan



Warehouse Operatives:
Jonathon & Gary

specula



Slightly oval medial end conforms to ear canal
Matt finish reduces glare
Serrated lateral end for easy grip

ear

			single use	re-usable		
Adept Specula 6mm	Pack of 100	NZ3230	£60.00	—	—	
Adept Specula 5mm	Pack of 100	NZ3231	£60.00	—	—	
Adept Specula 4mm	Pack of 100	NZ3232	£60.00	—	—	
Welch Allyn Child 2.75mm	Pack of 850	P60-100	£36.00	—	—	
Welch Allyn Adult 4.25mm	Pack of 850	P60-101	£36.00	—	—	
Gruber 1	Each	TE10-135	£2.50	E10-135	£14.00	
Gruber 2	Each	TE10-136	£2.50	E10-136	£14.00	
Gruber 3	Each	TE10-137	£2.50	E10-137	£14.00	
Gruber 4	Each	TE10-138	£2.50	E10-138	£14.00	
Tumarkin 1	Each	TE10-251	£2.50	E10-251	£25.00	
Tumarkin 2	Each	TE10-252	£2.50	E10-252	£25.00	
Tumarkin 3	Each	TE10-253	£2.50	E10-253	£25.00	
Tumarkin 4	Each	TE10-254	£2.50	E10-254	£25.00	
Tumarkin 5	Each	TE10-255	£2.50	E10-255	£25.00	
Tumarkin 6	Each	TE10-256	£2.50	E10-256	£25.00	
Tumarkin 7	Each	TE10-257	£2.50	E10-257	£25.00	
Tumarkin 8	Each	TE10-258	£2.50	E10-258	£25.00	



nasal

			single use	re-usable		
Thudicum 1	Each	TE12-124	£1.95	E12-124	£11.50	
Thudicum 2	Each	TE12-125	£1.95	E12-125	£11.50	
Thudicum 3	Each	TE12-126	£1.95	E12-126	£11.50	
Thudicum 4	Each	TE12-130	£1.95	E12-130	£11.50	
Thudicum 5	Each	TE12-131	£1.95	E12-131	£11.50	
Thudicum 6	Each	TE12-132	£1.95	E12-132	£11.50	
Thudicum 7	Each	TE12-133	£1.95	E12-133	£11.50	



Please contact our team on the following details:

Sheffmed Trade Services

Unit A, 40 Clifton Street
Sheffield, S9 2DQ
tel: 01142 617 161
fax: 01142 610 161
e-mail: sales@sheffmed.com

out of office hours contact:
Alan: 07711 324 645
Phill: 07947 128 347

suction tubes

Our unique range of suction tubes, fine ends and ointment syringes are manufactured to the highest standards. Packed sterile they combine safety with convenience

single use

P20-100	£55.00	Zoellner suction tube 30 degree no hole	Box of 50
P20-102	£59.00	Zoellner suction tube 30 degree with hole	Box of 50
P20-103	£55.00	Zoellner suction tube 30 degree no hole 45mm	Box of 50
P20-104	£55.00	Zoellner suction tube 30 degree with hole 45mm	Box of 50
P20-200	£39.00	Fine end 18 gauge	Box of 50
P20-201	£39.00	Fine end 20 gauge	Box of 50
P20-202	£39.00	Fine end 22 gauge	Box of 50
P20-120	£95.00	Frazier suction tube 3mm x 140mm	Box of 50
P20-121	£95.00	Frazier suction tube 2mm x 100mm	Box of 50
P20-122	£95.00	Frazier suction tube 3mm x 100mm	Box of 50
P20-510	£95.00	Magill suction tube 3mm fine 110m length	Box of 50
P20-125	£95.00	Bellucci suction tube 3mm fine 110mm length	Box of 50
P20-300	£147.00	Sinus fess malleable 2mm fine 110mm length	Box of 50
P20-400ST	£85.00	Long suction tube 3mm x 250mm	Box of 25
P20-410ST	£95.00	Long suction tube 3mm x 450mm	Box of 25
P20-420ST	£115.00	Long suction tube 3mm x 450mm	Box of 25
P20-430ST	£85.00	Long suction tube 4mm x 250mm	Box of 25
P20-440ST	£95.00	Long suction tube 4mm x 350mm	Box of 25
P20-450ST	£115.00	Long suction tube 4mm x 450mm	Box of 25

ointment tips

TE10-275S	£39.00	Ointment syringe with tip	Box of 50
TE10-276S	£29.00	Ointment tip only	Box of 50

metal suction tubes

		single use	re-usable
Metal fess malleable	Each	TE16-190S £3.95	E16-190 £35.00
Zoellners	Each	TG35-147S £2.90	G35-147 £17.25
Schucknecht 18g	Each	TG35-155S £3.90	G35-155 £22.50
Schucknecht 20g	Each	TG35-138S £3.90	G35-138 £22.50
Schucknecht 22g	Each	TG35-154S £3.90	G35-154 £22.50
Schucknecht 24g	Each	TG35-139S £3.90	G35-139 £22.50

instruments




			single use	re-usable
Crocodile forcep micro 3mm	Each	S11-330S	£7.90	E11-245 £120.00
Crocodile forcep standard 5mm	Each	TE11-241S	£7.90	E11-241 £120.00
Crocodile forcep large 7mm	Each	TE10-216S	£7.90	E10-216 £120.00
Henckel forcep	Each	S11-338S	£8.35	— —
Tilley forcep	Each	TE10-213S	£3.25	E10-213 £35.00
Micro scissors 5mm straight	Each	S11-345S	£11.90	E11-259 £120.00
Jobson horne metal	Each	TG24-135S	£1.30	G24-135 £12.75
Jobson horne carbon fibre	Pack of 100	P50-100	£45.00	— —
Myringotome straight (malleable)	Pack of 10	P70-101	£49.00	— —
Wax hook St Barts hook blunt	Pack of 100	NZ3441-100	£49.00	— —
Myringotomy set	Each	P70-100	£29.00	— —

Needle angled	Each	TE11-204S	£2.95	E11-204	£25.00
Needle curved	Each	TE11-202S	£2.95	E11-202	£25.00
Needle straight	Each	TE11-207S	£2.95	E11-207	£25.00
House curette	Each	TE11-215S	£3.95	E11-215	£25.00
Negus apirator	Each	TE12-183S	£3.95	E12-183	£25.00
Tracheal dilator adult	Each	TE19-102S	£4.95	E19-102	£35.00
Tracheal dilator child	Each	TE19-102CS	£4.95	E19-102C	£35.00
Tongue depressor 180mm	Each	TG18-108	£0.95	G18-108	£9.50
Lacks tongue depressor 13mm (Child)	Each	TG18-105	£2.00	G18-105	£9.75
Lacks tongue depressor 22mm (Adult)	Each	TG18-107	£2.00	G18-107	£10.50
Dental syringe 1.8cc	Each	TSP0068S	£5.95	SP0068	£44.50
Dental syringe 2.2cc	Each	TS11-229S	£6.95	S11-229	£46.25

single use re-usable

nasal splints

NZ3243	£98.00	Pack of 20
Malleable aluminium splint		
Can be trimmed to size		
Dorsal sponge support		
NZ3245	£40.00	Pack of 10
Anatomically designed to fit & support the nasal septum		
Pre-punched suture holes		
Saddle hinge with 8mm or 9mm gap		



headlights

micro spot headlight



High performance fibre optic headlight supplied with or without cross-band. Neoprene bolster provides excellent comfort for prolonged use.
Cable length 2.2m

08-370
08-375

Without cross-band **£950.00**
With cross-band **£950.00**

xenon headlight



08-137
08-137S



With integrated battery box. Very bright Xenon Headlight with high power optical system and fixed focus. 6000 candle power from 2 AA batteries

Rigid headband **£95.00**
Soft headband **£49.00**

classic headlights



Classic ENT headlight by Clar. 6 volt 2.4 watt with silver coated concave mirror.
Adjustable light focus
Comfortable headband with adjustment wheel
Weight 200g

08-13106 55mm **£150.00**
08-14106 98mm **£150.00**

binocular loupes



The high quality, coated optics ensure high resolution and colour authenticity. The angle of convergence of the light weight telescopic spectacles can be aligned independently of the inter-pupil distance and allows the optimal adjustment.

The flip-up magnifying unit enables spatial orientation during work pauses. The 3.5mm fibre bundle of F/O - illumination is guiding the light loss free to the light outlet. The high performance double lens system effect a homogenous spot which is adjustable in direction and size.

Both come with headband and fibre optic illumination
200-200 magnification 1.8x **£955.00**
200-205 magnification 2.5x **£955.00**

otoscopes



Fibre optic illumination with fitting for 2.5 and 4mm specula.

01-2060 Otoscope heine fit fibre optic **£70.00**

specula

01-9351 2.5mm pack of 500 **£95.00**
01-9361 4mm pack of 500 **£49.00**

headlight system complete

Halo150 comes complete with a Micro Spot headlight, 150w high performance halogen light source and mobile stand. (Xenon configuration also available)

12-750 **£1950.00**

